SEMESTER VI

INORGANIC CHEMISTRY

VI

Inorganic Preparations

Experiment No. 1 (3 & 6 Units)

Tris (Acetylacetonato) Iron (III), Fe(acac),

Aim: To prepare tris (acetylacetonato) iron (III).

Theory: This experiment is based on Green Chemistry Principles. It involves high atom economy, less hazardous chemical synthesis, environmentally clean and safe reaction at room temperature and minimises waste.

Metal acetyl acetonatos are highly efficient catalysts for a wide variety of organic transformations such as polymerisation, hydrogenation, coupling etc. They are used in rubber technology for vulcanization, for extraction and separation of metals, as fungicides, in pigments as colour stabilizers, as carbon scavengers for diesel fuels, as combustion control catalysts for rocket fuels etc.

Tris (acetyl acetonato) iron (III) is prepared directly from Fe(OH)₃ and acetyl acetone reagent without the use of any buffer like sodium acetate which contaminates the product.

Requirements: FeCl₃ • 6H₂O, 1 N KOH solution, acetyl acetone reagent, buchner funnel, ice bath, fused CaCl₂, Whatman filter paper No. 1, vacuum desiccator, distilled water, 100 cm³ beaker, etc.

Reaction: $Fe(OH)_3 + 3CH_3COCH_2COCH_3 \rightarrow Fe(CH_2COCH_2COCH_3)_3 + 3H_2O$

Procedure:

- (1) Dissolve W g (i.e. 1.5 g or the supplied quantity) of FeCl₃•6H₂O in 10 cm³ of distilled water in a 100 cm³ beaker. (Heat it if necessary).
- (2) Add about 30 cm³ of 1N KOH solution in small volumes to the above FeCl₃ solution. Stir the solution after each addition of KOH solution to precipitate Fe(OH)₃ completely.
- (3) Allow the precipitate to settle down completely. Supernatant solution becomes colourless
- (4) Decant the supernatant colourless solution and wash the precipitate 2-3 times with distilled water using 5-7 cm³ each time. Decant the solution each time after washing.
- (5) Filter the precipitate of Fe(OH)₃ through buchner fuunel containing Whatman filter No. 1, wash it twice with distilled water.
- (6) Transfer the precipitate to 100 cm³ beaker and keep it in an icebath for 10 minutes.
- (7) Add about 4-5 cm³ of distilled acetylacetone reagent to the cold Fe(OH)₃ precipitate and stir it well with a glass rod.
- (8) Allow the reaction mixture to stand at room temperature for about 30 minutes. (Stir it occassionally with the glass rod).
- (9) Shiny deep red coloured iron-acetyl acetonato complex gradually forms. Keep the beaker in an icebath for about 15 minutes.
- (10) Filter the complex through buchner funnel containing Whatman filter paper No. 41 and dry it in vacuum over fused CaCl₂ in a desicçator.
- (11) Weigh the complex to obtain the yield (x g).

Observations and Calculations:

Weight of the complex $Fe(acac)_3 = x g$.

Yield of the complex:

- (A) Observed Yield of the complex = x g.
- (B) Theoretical Yield of the complex:

Now, 270.35 g of FeCl₃•6H₂O gives 355.85 g of Fe(acac)₃

Experiment No. 2 (6 Units Only)

137